

STOCK IMMERSION STERILIZERS

1300mm MRI MULTI-PROCESS

- Process Modes: full immersion, partial immersion, spray
- Agitation Modes: static, rotation (E-O-E), oscillation
- Process/storage vessels, piping, panel and rotor constructed of mild or 304L stainless steel (316 SS optional)
- Insulated shell with mirror polished SS cladding
- Easily integrates with automated shuttles for loading and unloading
- Includes manual door, optional automatic vertical door and individual pneumatic cage clampdown
- Full and partial-immersion includes direct heating and cooling. (Indirect heating and cooling optional)
- Forced convection is obtained with rotary operation accelerating product heat transfer to the coldest spot in the containers, even with varying viscosity products
- Short process times preserve natural quality, vitamins and nutritional values
- Agitation prevents burn and provides uniform processes
- Optimum temperature distribution (+/- 1°F) and accurate pressure control (+/- 1 psig)

